

Implementation of software and hardware of free access for the technological solution to illustrate physical phenomena in the classroom in a didactically and pedagogical way

For more information and laboratory setups.


Contact:	Diego Parra	Lic. in physics,
	Luis Ali Ortiz	P.C.I,
	Juan Bastidas	P.C.I,
	Miguel Gonzalez	Ing. systems

Email:	dparra@opensai.org
Email:	lortiz@opensai.org
Email:	jbastidas@opensai.org
Email:	mgonzales@opensai.org

Vehicle

This vehicle is hardware to free access and this is controlled for android and GNU-Linux platform or smart glass BT-200 moveiro device to EPSON.

This software is to free access and development in a python language.




Ultra-sound sensor

data acquisition card
with radio frequency
and bluetooth
communication

Sensor based to
passive electro-
reception

VELCRO

Sensor based to
Active electro-
reception




APP

- FREinfraROSI is software to free access, writing this app in a python language and build this in a GNU-Linux platform.
- With this app, can development physical laboratory, but example free harmonic movement, calculate to decay coefficient to the flux, Law of decay radiation flux density, and visualize diffraction experiment.

• free harmonic movement

Sensor and move car button

Exit button



Motorize vehicle

Law of decay
radiation flux
density

Diffraction button

Capture image
button

Clear Button

Sensor and move car button



measure
length button

Pasive electro-
reception
button

Light
sensor
button

Move car to
righ button

Active electro-
reception
button

Move car to
righ button

Stop button

Finish button

length of one rope from a count the oscillations

Put data acquisition card with sensor based to active electro-reception and the system mass-rope in front, use a software FreePS in a device android or a smart glass BT-200 EPSON.

FreePS software begin to communication bluetooth with a acquisition card, and count a oscillation to system mass in a rope.

$$(T^2) / (4 * \pi^2) = (L / g);$$

$$g / (4 * f^2 * \pi^2) = L$$


g = gravity acceleration

$\pi = 3.14159265$

T = oscillation period

L = length string in meter

f = Frequency



Law of decay radiation flux density

Put a mirror in front to the vehicle at a distance of 8.5 cm, the vehicle move 2 mm parallel to the radiation flux, which flux decays with the inverse distance squared to distance of separation to the source. Acquisition card with the sensor active electro-reception, capture the data and sends it to the device, which uses a logarithmic linear regression to find the exponent of radiation decay because of changing in the distance and calculates a density flux of radiation means emitted by the diode.

Calculate exponent and flux density

$$E = d\Phi/ds$$

$$I = d\Phi/d\omega$$

$$d\omega = ds/r^2$$

$$E = I/r^2$$


Statistics

$$\ln(Y^*) = \ln(a) + b\ln(X) \cdots \implies V^* = A + bU,$$

$$\text{Exponent} \quad b = \frac{S_{UV}}{S_U^2} = \frac{\frac{1}{n} \sum_{i=1}^n UV - \bar{U}\bar{V}}{\frac{1}{n} \sum_{i=1}^n U^2 - \bar{U}^2},$$

$$\text{Statistics last} \quad Y^* = aX^b = [\text{antiln}(\bar{V} - b\bar{U})] * X^{\left(\frac{\frac{1}{n} \sum_{i=1}^n UV - \bar{U}\bar{V}}{\frac{1}{n} \sum_{i=1}^n U^2 - \bar{U}^2}\right)},$$

Software



THANKS

AGRADECIMIENTOS

Agradeciendo al grupo de investigación en física e informática 'FISINFOR' de la Universidad Distrital Francisco José de Caldas, al equipo de HACK-VISION, al equipo de OpenSAI, a los usuarios de este software y hardware, y de manera cariñosa a todas las personas que contribuyen para que el conocimiento sea libre.

Agradecimientos especiales:

* PhD. Julián Andrés Salamanca Bernal